

IVUS Plus

Complement your ClearView™ with Longitudinal and Digital Review

Major Benefits:

- ◆ Adds PC technology to **modernize your existing IVUS system**
- ◆ Adds **digital review and storage**
- no more hunting back and forth on the video tape
- ◆ **Real-time rotation of longitudinal view** for easier identification of region of interest and lesion length
- ◆ **DICOM Read and Write (Optional)**
- ◆ **2D and 3D volumetric analysis**
- ◆ **Reporting and networking capabilities**
- ◆ **Automatic export** of numerical results to spreadsheet file
- ◆ **Make AVI files** of all or part of a pullback



Update your existing Boston Scientific ClearView™ IVUS system with INDEC System's IVUS Plus. The PC-based computer system, with Windows 2000 operating system and proven cardiovascular software solutions, modernizes and totally updates your current IVUS capabilities with the latest in clinical and research features. The IVUS Plus upgrade provides real time longitudinal view, digital storage and review, volume measurements, true mutli-tasking and network connectivity. No more bulky video tape required. Review and store, process and analyze your IVUS pullbacks to disk, CD ROM media or a network. In addition, INDEC's IVUS Plus provides the ability to browse the web, generate reports, read and write DICOM files, make AVI files and easily interface to PowerPoint™.

IVUS Plus provides the enhanced electronics and proven IVUS applications software to modernize your ClearView system.



Enhance your ClearView system by adding PC technology and easy access for USB, IEEE1394, printer and networking connectivity.

IVUS Plus adds a modern flat screen monitor which can be easily folded out of the way when not in use.



Longitudinal view is acquired in real-time so it is immediately available during the procedure.

Galaxy Ready: when your ready for a new IVUS system, this upgrade can be easily reversed and you then have a stand alone analysis workstation that can work with your old files on CD or video tape or files from your new Galaxy machine. Installation and Training take place at your hospital site - no need to send your system anywhere to be upgraded. Training also includes a "Quick Start" video library of how to use the key features. Please visit our web site - www.indecsystems.com for further information.

510K Clearance Received

INDEC Systems, Inc.

Tel: 1-650-903-9755

Fax: 1-650-919-0150

www.indecsystems.com

IVPI-0820031-C