



Maximum  
Exposure to Make  
Valve Procedures  
Easier



ESTECH'S Valve Attachments work in harmony with the Universal Stabilizer platform to provide optimal surgical exposure for valve procedures.

**VALVE**

# VALVE

## Atrial Rakes

ESTECH'S Atrial Rakes coupled with the Universal Stabilizer platform provide optimal surgical exposure for mitral valve procedures.



**Atrial Rake**  
Small



**Atrial Rake**  
Medium



**Atrial Rake**  
Large

## Aortic Valve Assistants — Craver Design

ESTECH'S Aortic Valve Assistants coupled with the Universal Stabilizer platform provide optimal surgical exposure for aortic valve procedures (also useful for mitral valve exposure when using a transeptal approach).




**Aortic Valve Assistant**  
Small



**Aortic Valve Assistant**  
Medium



**Aortic Valve Assistant**  
Large



## Universal Stabilizer Arm HP™

Low-profile articulating arm mounts to sternal retractor, providing superior positional ability and strength.

Product Description	Catalog Number
Universal Stabilizer Arm HP	401-152HP
Atrial Rake— Small	400-441
Atrial Rake— Medium	400-442
Atrial Rake— Large	400-443
Aortic Valve Assistant— Small	400-451
Aortic Valve Assistant— Medium	400-452
Aortic Valve Assistant— Large	400-453

11566-LIT-A  
Patent No. 6,860,668  
Other Patents Pending

