


Cardiac catheterization is performed to diagnose the nature, severity, and precise location of blockages and other abnormalities of the heart. This technique represents the most essential diagnostic tool in the management of patients with cardiovascular disease.




## Diagnostic Cardiology Catheters


DIAGNOSTIC CARDIOLOGY CATHETERS

**MERTMEDICA**


- a) JUDKINS LEFT 3.5
- b) JUDKINS LEFT 4.0
- c) JUDKINS LEFT 5.0
- d) JUDKINS LEFT 6.0
- e) JUDKINS LEFT SHORT TIP 3.5
- f) JUDKINS LEFT SHORT TIP 4.0
- g) JUDKINS LEFT SHORT TIP 5.0
- h) JUDKINS LEFT SHORT TIP 6.0




- a) AMPLATZ LEFT 1
- b) AMPLATZ LEFT 2
- c) AMPLATZ LEFT 3
- d) AMPLATZ LEFT 4
- e) AMPLATZ RIGHT 1
- f) AMPLATZ RIGHT 2
- g) AMPLATZ RIGHT 3
- h) AMPLATZ RIGHT MODIFIED
- i) PROGRESSIVE RIGHT (non torque)




- a) JUDKINS RIGHT 3.5
- b) JUDKINS RIGHT 4.0
- c) JUDKINS RIGHT 5.0
- d) JUDKINS RIGHT 6.0
- e) JUDKINS RIGHT SHORT TIP 3.5
- f) JUDKINS RIGHT SHORT TIP 4.0
- g) JUDKINS RIGHT SHORT TIP 5.0
- h) JUDKINS RIGHT SHORT TIP 6.0
- i) JUDKINS RIGHT MODIFIED 3.5
- j) JUDKINS RIGHT MODIFIED 4.0
- k) JUDKINS RIGHT MODIFIED 5.0
- l) JUDKINS RIGHT MODIFIED 6.0




- a) STRAIGHT PIGTAIL
- b) 145° PIGTAIL
- c) 155° PIGTAIL
- d) MULTIPURPOSE A1
- e) MULTIPURPOSE A2
- f) MULTIPURPOSE B1
- g) MULTIPURPOSE B2
- h) SONES




- a) LEFT CORONARY BYPASS
- b) RIGHT CORONARY BYPASS
- c) INTERNAL MAMMARY



- a) INTERNAL MAMMARY  
LEFT BRACHIAL/RADIAL 1
- b) INTERNAL MAMMARY  
LEFT BRACHIAL/RADIAL 2
- c) LEFT CORONARY  
BRACHIAL/RADIAL 1
- d) LEFT CORONARY  
BRACHIAL/RADIAL 2
- e) LEFT CORONARY  
BRACHIAL/RADIAL 3
- f) RIGHT CORONARY  
BRACHIAL/RADIAL
- g) AMPLATZ LEFT CORONARY  
BRACHIAL/RADIAL



- a) PED. JUDKINS LEFT 1.5
- b) PED. JUDKINS LEFT 2.0
- c) PED. JUDKINS LEFT 2.5
- d) PED. JUDKINS LEFT 3.0
- e) PED. JUDKINS RIGHT 1.5
- f) PED. JUDKINS RIGHT 2.0
- g) PED. JUDKINS RIGHT 2.5
- h) PED. JUDKINS RIGHT 3.0





**DIAGNOSTIC CARDIOLOGY CATHETERS B**


**PERFORMA® & SOFTOUCH®**

Judkins Left	B.1 & 3
Judkins Left Short Tip	B.2, 3
Judkins Right	B.4 & 6
Judkins Right Short Tip	B.5, 6
Judkins Right Modified	B.5, 6
Amplatz Left	B.7
Amplatz Right	B.8
Left Coronary Bypass	B.8
Right Coronary Bypass	B.9
Progressive Right	B.9
Internal Mammary	B.9
Cardiac Pigtails	B.10, 11, 12
Pediatric / Infant Pigtails	B.13
Pediatric Cardiology	B.14
Multipurpose	B.15, 16
Sones	B.16
Brachial / Radial	B.17
Cardiac MultiPack / Cardiac MultiPack Plus	B.18, 19

catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--

### PERFORMA® JUDKINS LEFT CATHETERS

<b>7700-10</b>	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL3.5	bumper tip
<b>7700-20</b>	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL4.0	bumper tip
<b>7700-30</b>	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL5.0	bumper tip
<b>7700-40</b>	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	JL6.0	bumper tip
<b>7521-13</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL3.5	
<b>7521-21</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL4.0	
<b>7521-33</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL5.0	
<b>7521-43</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL6.0	
<b>7521-F1</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JL4.5	
<b>7521-53</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL3.5	
<b>7521-63</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL4.0	
<b>7521-73</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL5.0	
<b>7521-83</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JL6.0	
<b>7601-13</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL3.5-High Flow	non-braided
<b>7601-21</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL4.0-High Flow	non-braided
<b>7601-33</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL5.0-High Flow	non-braided
<b>7601-43</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL6.0-High Flow	non-braided
<b>7501-13</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL3.5	
<b>7501-21</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL4.0	
<b>7501-H1</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL4.5	
<b>7501-33</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL5.0	
<b>7501-43</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL6.0	
<b>7501-53</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL3.5	
<b>7501-63</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL4.0	
<b>7501-73</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL5.0	
<b>7501-83</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JL6.0	
<b>7574-13</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL3.5-High Flow	non-braided
<b>7574-21</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL4.0-High Flow	non-braided
<b>7574-33</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL5.0-High Flow	non-braided
<b>7574-43</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL6.0-High Flow	non-braided
<b>7670-13</b>	7	.064"/1.62mm	100cm	7.2cm	.038" (0.97mm)	JL3.5	
<b>7670-21</b>	7	.064"/1.62mm	100cm	7.2cm	.038" (0.97mm)	JL4.0	
<b>7670-33</b>	7	.064"/1.62mm	100cm	8.0cm	.038" (0.97mm)	JL5.0	
<b>7670-43</b>	7	.064"/1.62mm	100cm	9.0cm	.038" (0.97mm)	JL6.0	



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description
----------------	-------------	-------------------	--------	---------------	------------------------	-------------

**PERFORMA® JUDKINS LEFT SHORT-TIP CATHETERS**

<b>7501-A3</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL3.5 ST
<b>7501-B1</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL4.0 ST
<b>7501-C3</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL5.0 ST
<b>7501-D3</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JL6.0 ST




DIAGNOSTIC CARDIOLOGY CATHETERS

catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description
----------------	-------------	-------------------	--------	---------------	------------------------	-------------

### SOFTOUCH® JUDKINS LEFT CATHETERS

<b>7402-13</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JL3.5	
<b>7402-21</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JL4.0	
<b>7402-33</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JL5.0	
<b>7402-43</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JL6.0	
<b>7431-13</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL3.5-High Flow	non-braided
<b>7431-21</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL4.0-High Flow	non-braided
<b>7431-33</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL5.0-High Flow	non-braided
<b>7431-43</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	JL6.0-High Flow	non-braided
<b>7402-53</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL3.5	
<b>7402-61</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL4.0	
<b>7402-73</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL5.0	
<b>7402-83</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL6.0	
<b>7221-13</b>	6	.051"/1.29mm	125cm	2.5cm	.038" (0.97mm)	JL3.5	
<b>7221-23</b>	6	.051"/1.29mm	125cm	2.5cm	.038" (0.97mm)	JL4.0	
<b>7221-33</b>	6	.051"/1.29mm	125cm	2.5cm	.038" (0.97mm)	JL5.0	
<b>7221-43</b>	6	.051"/1.29mm	125cm	2.5cm	.038" (0.97mm)	JL6.0	
<b>7431-53</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL 3.5-High Flow	non-braided
<b>7431-61</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL 4.0-High Flow	non-braided
<b>7431-73</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL 5.0-High Flow	non-braided
<b>7431-83</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	JL 6.0-High Flow	non-braided



### SOFTOUCH® JUDKINS LEFT SHORT-TIP CATHETERS


<b>7413-13</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL3.5 ST
<b>7413-21</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL4.0 ST
<b>7413-33</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL5.0 ST
<b>7413-43</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JL6.0 ST



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description
----------------	-------------	-------------------	--------	---------------	------------------------	-------------

### PERFORMA® JUDKINS RIGHT CATHETERS

<b>7701-10</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR3.5 bumper tip
<b>7701-20</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR4.0 bumper tip
<b>7701-30</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR5.0 bumper tip
<b>7701-40</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	JR6.0 bumper tip
<b>7523-13</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR3.5
<b>7523-21</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR4.0
<b>7523-33</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR5.0
<b>7523-43</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR6.0
<b>7523-53</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR3.5
<b>7523-63</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR4.0
<b>7523-73</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR5.0
<b>7523-83</b>	5	.046"/1.17mm	125cm	2.5cm	.038" (0.97mm)	JR6.0
<b>7503-13</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR3.5
<b>7503-21</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR4.0
<b>7503-33</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR5.0
<b>7503-43</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR6.0
<b>7503-53</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR3.5
<b>7503-63</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR4.0
<b>7503-73</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR5.0
<b>7503-83</b>	6	.054"/1.37mm	125cm	1.9cm	.038" (0.97mm)	JR6.0
<b>7661-13</b>	7	.064"/1.62mm	100cm	2.4cm	.038" (0.97mm)	JR3.5
<b>7661-21</b>	7	.064"/1.62mm	100cm	2.1cm	.038" (0.97mm)	JR4.0
<b>7661-33</b>	7	.064"/1.62mm	100cm	2.4cm	.038" (0.97mm)	JR5.0
<b>7661-43</b>	7	.064"/1.62mm	100cm	2.4cm	.038" (0.97mm)	JR6.0
<b>7661-53</b>	7	.064"/1.62mm	125cm	2.4cm	.038" (0.97mm)	JR3.5
<b>7661-63</b>	7	.064"/1.62mm	125cm	2.1cm	.038" (0.97mm)	JR4.0
<b>7661-73</b>	7	.064"/1.62mm	125cm	2.4cm	.038" (0.97mm)	JR5.0
<b>7661-83</b>	7	.064"/1.62mm	125cm	2.4cm	.038" (0.97mm)	JR6.0



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description
----------------	-------------	-------------------	--------	---------------	------------------------	-------------

**PERFORMA® JUDKINS RIGHT SHORT-TIP CATHETERS**

<b>7523-E3</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR3.5 ST
<b>7523-F1</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR4.0 ST
<b>7523-G3</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR5.0 ST
<b>7523-H3</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	JR6.0 ST
<b>7503-A3</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR3.5 ST
<b>7503-B1</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR4.0 ST
<b>7503-C3</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR5.0 ST
<b>7503-D3</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR6.0 ST



**PERFORMA® JUDKINS RIGHT MODIFIED CATHETERS**


<b>7502-13</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR3.5 MOD
<b>7502-23</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR4.0 MOD
<b>7502-33</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR5.0 MOD
<b>7502-43</b>	6	.054"/1.37mm	100cm	1.9cm	.038" (0.97mm)	JR6.0 MOD



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description
----------------	-------------	-------------------	--------	---------------	------------------------	-------------

**SOFTOUCH® JUDKINS RIGHT CATHETERS**

<b>7401-P3</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR3.5
<b>7401-Q1</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR4.0
<b>7401-R3</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR5.0
<b>7401-S3</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR6.0
<b>7401-A3</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR3.5
<b>7401-B1</b>	6	.051"/1.29mm	100cm	2.4cm	.038" (0.97mm)	JR4.0
<b>7401-C3</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR5.0
<b>7401-D3</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR6.0
<b>7220-13</b>	6	.051"/1.29mm	125cm	2.5cm	.038" (0.97mm)	JR3.5
<b>7220-23</b>	6	.051"/1.29mm	125cm	2.4cm	.038" (0.97mm)	JR4.0
<b>7220-33</b>	6	.051"/1.29mm	125cm	2.4cm	.038" (0.97mm)	JR5.0
<b>7220-43</b>	6	.051"/1.29mm	125cm	2.4cm	.038" (0.97mm)	JR6.0




**SOFTOUCH® JUDKINS RIGHT SHORT-TIP CATHETERS**

<b>7411-13</b>	6	.051"/1.29mm	100cm	1.9cm	.038" (0.97mm)	JR3.5 ST
<b>7411-21</b>	6	.051"/1.29mm	100cm	1.9cm	.038" (0.97mm)	JR4.0 ST
<b>7411-33</b>	6	.051"/1.29mm	100cm	1.9cm	.038" (0.97mm)	JR5.0 ST
<b>7411-43</b>	6	.051"/1.29mm	100cm	1.9cm	.038" (0.97mm)	JR6.0 ST



**SOFTOUCH® JUDKINS RIGHT MODIFIED CATHETERS**

<b>7401-13</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR3.5 MOD
<b>7401-21</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR4.0 MOD
<b>7401-33</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR5.0 MOD
<b>7401-43</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	JR6.0 MOD
<b>7401-53</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR3.5 MOD
<b>7401-61</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR4.0 MOD
<b>7401-73</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR5.0 MOD
<b>7401-83</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	JR6.0 MOD




catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description
----------------	-------------	-------------------	--------	---------------	------------------------	-------------

### PERFORMA® AMPLATZ LEFT CATHETERS

<b>7703-10</b>	4	.042"/1.07mm	100cm	5.0cm	.038" (0.97mm)	AL1 bumper tip
<b>7703-20</b>	4	.042"/1.07mm	100cm	2.2cm	.038" (0.97mm)	AL2 bumper tip
<b>7703-30</b>	4	.042"/1.07mm	100cm	2.3cm	.038" (0.97mm)	AL3 bumper tip
<b>7703-40</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	AL4 bumper tip
<b>7527-13</b>	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AL1
<b>7527-23</b>	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AL2
<b>7527-33</b>	5	.046"/1.17mm	100cm	2.0cm	.038" (0.97mm)	AL3
<b>7527-43</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	AL4
<b>7507-13</b>	6	.054"/1.37mm	100cm	5.0cm	.038" (0.97mm)	AL1
<b>7507-23</b>	6	.054"/1.37mm	100cm	5.0cm	.038" (0.97mm)	AL2
<b>7507-33</b>	6	.054"/1.37mm	100cm	2.4cm	.038" (0.97mm)	AL3
<b>7507-43</b>	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	AL4
<b>7662-13</b>	7	.064"/1.62mm	100cm	6.5cm	.038" (0.97mm)	AL1
<b>7662-23</b>	7	.064"/1.62mm	100cm	7.5cm	.038" (0.97mm)	AL2
<b>7662-33</b>	7	.064"/1.62mm	100cm	7.5cm	.038" (0.97mm)	AL3
<b>7662-43</b>	7	.064"/1.62mm	100cm	7.5cm	.038" (0.97mm)	AL4

### SOFTOUCH® AMPLATZ LEFT CATHETERS

<b>7406-13</b>	5	.044"/1.12mm	100cm	5.0cm	.038" (0.97mm)	AL1
<b>7406-23</b>	5	.044"/1.12mm	100cm	5.0cm	.038" (0.97mm)	AL2
<b>7406-33</b>	5	.044"/1.12mm	100cm	2.0cm	.038" (0.97mm)	AL3
<b>7406-43</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	AL4
<b>7406-53</b>	6	.051"/1.29mm	100cm	5.0cm	.038" (0.97mm)	AL1
<b>7406-63</b>	6	.051"/1.29mm	100cm	5.0cm	.038" (0.97mm)	AL2
<b>7406-73</b>	6	.051"/1.29mm	100cm	2.4cm	.038" (0.97mm)	AL3
<b>7406-83</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	AL4




catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--

**PERFORMA® AMPLATZ RIGHT CATHETERS**

<b>7702-40</b>	4	.042"/1.07mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	bumper tip
<b>7702-10</b>	4	.042"/1.07mm	100cm	3.3cm	.038" (0.97mm)	AR1	bumper tip
<b>7702-20</b>	4	.042"/1.07mm	100cm	5.0cm	.038" (0.97mm)	AR2	bumper tip
<b>7702-30</b>	4	.042"/1.07mm	100cm	5.0cm	.038" (0.97mm)	AR3	bumper tip
<b>7526-11</b>	5	.046"/1.17mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	
<b>7528-13</b>	5	.046"/1.17mm	100cm	3.3cm	.038" (0.97mm)	AR1	
<b>7528-23</b>	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AR2	
<b>7528-33</b>	5	.046"/1.17mm	100cm	5.0cm	.038" (0.97mm)	AR3	
<b>7506-11</b>	6	.054"/1.37mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	
<b>7508-13</b>	6	.054"/1.37mm	100cm	3.0cm	.038" (0.97mm)	AR1	
<b>7508-23</b>	6	.054"/1.37mm	100cm	3.5cm	.038" (0.97mm)	AR2	
<b>7508-33</b>	6	.054"/1.37mm	100cm	4.0cm	.038" (0.97mm)	AR3	
<b>7663-41</b>	7	.064"/1.62mm	100cm	4.0cm	.038" (0.97mm)	AR MOD	
<b>7663-13</b>	7	.064"/1.62mm	100cm	4.0cm	.038" (0.97mm)	AR1	
<b>7663-23</b>	7	.064"/1.62mm	100cm	5.0cm	.038" (0.97mm)	AR2	
<b>7663-33</b>	7	.064"/1.62mm	100cm	6.0cm	.038" (0.97mm)	AR3	

**SOFTOUCH® AMPLATZ RIGHT CATHETERS**

<b>7404-11</b>	5	.044"/1.12mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	
<b>7404-21</b>	6	.051"/1.29mm	100cm	3.3cm	.038" (0.97mm)	AR MOD	
<b>7403-43</b>	6	.051"/1.29mm	100cm	3.0cm	.038" (0.97mm)	AR1	
<b>7403-53</b>	6	.051"/1.29mm	100cm	3.5cm	.038" (0.97mm)	AR2	
<b>7403-63</b>	6	.051"/1.29mm	100cm	4.0cm	.038" (0.97mm)	AR3	



**PERFORMA® LEFT CORONARY BYPASS CATHETERS**

<b>7709-10</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	LCB	bumper tip
<b>7534-13</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	LCB	
<b>7510-13</b>	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	LCB	
<b>7665-13</b>	7	.064"/1.62mm	100cm	9.0cm	.038" (0.97mm)	LCB	

**SOFTOUCH® LEFT CORONARY BYPASS CATHETERS**

<b>7851-23</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	LCB	
<b>7399-23</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	LCB	



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--

### PERFORMA® RIGHT CORONARY BYPASS CATHETERS

<b>7709-20</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	RCB	bumper tip
<b>7534-23</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	RCB	
<b>7510-23</b>	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	RCB	
<b>7665-23</b>	7	.064"/1.62mm	100cm	7.5cm	.038" (0.97mm)	RCB	

### SOFTOUCH® RIGHT CORONARY BYPASS CATHETERS

<b>7851-13</b>	5	.044"/1.12mm	100cm	2.5cm	.038" (0.97mm)	RCB	
<b>7399-13</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	RCB	



### PERFORMA® PROGRESSIVE RIGHT CORONARY CATHETERS (NON-TORQUE)

<b>7707-10</b>	4	.042"/1.07mm	100cm	2.5cm	.038" (0.97mm)	PRC	bumper tip
<b>7600-11</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	PRC	non-braided
<b>7532-13</b>	5	.046"/1.17mm	100cm	2.5cm	.038" (0.97mm)	PRC	
<b>7572-11</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	PRC	non-braided
<b>7512-13</b>	6	.054"/1.37mm	100cm	2.5cm	.038" (0.97mm)	PRC	
<b>7512-23</b>	6	.059"/1.49mm	125cm	2.5cm	.038" (0.97mm)	PRC	
<b>7668-13</b>	7	.064"/1.62mm	100cm	2.5cm	.038" (0.97mm)	PRC	

### SOFTOUCH® PROGRESSIVE RIGHT CORONARY CATHETERS (NON-TORQUE)

<b>7430-11</b>	5	.052"/1.32mm	100cm	2.5cm	.038" (0.97mm)	PRC	non-braided
<b>7430-21</b>	6	.059"/1.49mm	100cm	2.5cm	.038" (0.97mm)	PRC	non-braided
<b>7488-13</b>	6	.051"/1.29mm	100cm	2.5cm	.038" (0.97mm)	PRC	
<b>7488-23</b>	6	.051"/1.29mm	125cm	2.5cm	.038" (0.97mm)	PRC	



### PERFORMA® INTERNAL MAMMARY CATHETERS

<b>7706-10</b>	4	.042"/1.07mm	100cm	3.3cm	.038" (0.97mm)	IMA	bumper tip
<b>7531-30</b>	5	.046"/1.17mm	65cm	3.3cm	.038" (0.97mm)	IMA	
<b>7531-13</b>	5	.046"/1.17mm	100cm	3.3cm	.038" (0.97mm)	IMA	
<b>7511-13</b>	6	.054"/1.37mm	100cm	3.3cm	.038" (0.97mm)	IMA	
<b>7666-13</b>	7	.064"/1.62mm	100cm	7.5cm	.038" (0.97mm)	IMA	

### SOFTOUCH® INTERNAL MAMMARY CATHETERS

<b>7852-13</b>	5	.044"/1.12mm	100cm	3.3cm	.038" (0.97mm)	IMA	
<b>7395-13</b>	6	.051"/1.29mm	100cm	3.3cm	.038" (0.97mm)	IMA	



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	#side ports	small pigtail curve
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	-------------	---------------------

**PERFORMA® STRAIGHT CARDIAC PIGTAIL CATHETERS**

<b>7704-70</b>	4	.042"/1.07mm	80 cm	6.5 cm	.038" (0.97mm)	Braided PIG	5	
<b>7704-10</b>	4	.042"/1.07mm	110 cm	6.5 cm	.038" (0.97mm)	Braided PIG	5	
<b>7776-11</b>	5	.052"/1.32mm	90cm	5.0cm	.038" (0.97mm)	Braided PIG	5	yes
<b>7776-21</b>	5	.052"/1.32mm	110cm	5.0cm	.038" (0.97mm)	Braided PIG	5	yes
<b>7776-B1</b>	5	.052"/1.32mm	110cm	4.6cm	.038" (0.97mm)	Non-Braided PIG	5	yes
<b>7602-11</b>	5	.052"/1.32mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided PIG	5	
<b>7529-11</b>	5	.046"/1.17mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7529-43</b>	5	.046"/1.17mm	125cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7569-11</b>	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided PIG	5	
<b>7509-11</b>	6	.054"/1.37mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7777-B1</b>	6	.059"/1.49mm	110cm	4.8cm	.038" (0.97mm)	Non-Braided PIG	5	yes
<b>7777-21</b>	6	.059"/1.49mm	110cm	4.6cm	.038" (0.97mm)	Braided PIG	5	yes
<b>7509-43</b>	6	.059"/1.49mm	125cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7569-43</b>	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided PIG	5	
<b>7669-11</b>	7	.064"/1.62mm	110cm	6.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7669-43</b>	7	.064"/1.62mm	125cm	6.5cm	.038" (0.97mm)	Braided PIG	5	



**PIG**

**SOFTOUCH® STRAIGHT CARDIAC PIGTAIL CATHETERS**

<b>7779-51</b>	5	.052"/1.32mm	110cm	4.7cm	.038" (0.97mm)	Non-Braided PIG	5	
<b>7779-G1</b>	5	.044"/1.12mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7779-A1</b>	6	.051"/1.29mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7779-P1</b>	6	.051"/1.29mm	110cm	8.0cm	.038" (0.97mm)	Mod. Braided PIG	10	
<b>7779-S1</b>	6	.051"/1.29mm	110cm	7.5cm	.038" (0.97mm)	Braided PIG	2/3*	
<b>7779-L1</b>	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided PIG	5	
<b>7780-13</b>	6	.051"/1.29mm	125cm	7.5cm	.038" (0.97mm)	Braided PIG	5	
<b>7222-13</b>	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided PIG	5	



**PIG**

\* 2 regular sized side ports and 3 larger sized side ports.

catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	#side ports	small pigtail curve	angle length
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	-------------	---------------------	--------------

### SOFTOUCH® CARDIAC PIGTAIL CATHETERS with lateral opposed sideports

<b>7223-11</b>	6	.051"/1.29mm	110cm	6.0cm	.038" (0.97mm)	Braided PIG	6		
<b>7223-23</b>	6	.051"/1.29mm	110cm	6.0cm	.038" (0.97mm)	Braided 145 PIG	6		
<b>7223-33</b>	6	.051"/1.29mm	110cm	6.0cm	.038" (0.97mm)	Braided 155 PIG	6		



### PERFORMA® 145 DEGREE CARDIAC PIGTAIL CATHETERS

<b>7704-20</b>	4	.042"/1.07mm	110cm	6.5cm	.038" (0.97mm)	Braided 145 PIG	5		5.5cm
<b>7776-C3</b>	5	.046"/1.17mm	110cm	4.6cm	.038" (0.97mm)	Non-Braided 145 PIG	5	yes	5.5cm
<b>7776-33</b>	5	.046"/1.17mm	110cm	5.0cm	.038" (0.97mm)	Braided 145 PIG	5	yes	5.5cm
<b>7529-53</b>	5	.046"/1.17mm	125cm	7.5cm	.038" (0.97mm)	Braided 145 PIG	5		7.5cm
<b>7602-23</b>	5	.052"/1.32mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided 145 PIG	5		7.5cm
<b>7529-23</b>	5	.046"/1.17mm	110cm	7.5cm	.038" (0.97mm)	Braided 145 PIG	5		7.5cm
<b>7569-23</b>	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided 145 PIG	5		5.5cm
<b>7509-23</b>	6	.054"/1.37mm	110cm	7.5cm	.038" (0.97mm)	Braided 145 PIG	5		5.5cm
<b>7777-C3</b>	6	.059"/1.49mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided 145 PIG	5	yes	5.5cm
<b>7777-33</b>	6	.054"/1.37mm	110cm	4.6cm	.038" (0.97mm)	Braided 145 PIG	5	yes	5.5cm
<b>7509-53</b>	6	.054"/1.37mm	125cm	7.5cm	.038" (0.97mm)	Braided 145 PIG	5		5.5cm
<b>7569-53</b>	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided 145 PIG	5		5.5cm
<b>7669-23</b>	7	.064"/1.62mm	110cm	6.5cm	.038" (0.97mm)	Braided 145 PIG	5		7.5cm
<b>7669-53</b>	7	.064"/1.62mm	125cm	6.5cm	.038" (0.97mm)	Braided 145 PIG	5		7.5cm

### SOFTOUCH® 145 DEGREE CARDIAC PIGTAIL CATHETERS

<b>7779-63</b>	5	.052"/1.32 mm	110cm	4.7cm	.038" (0.97mm)	Non-Braided 145 PIG	5		5.5cm
<b>7779-B3</b>	6	.051"/1.29mm	110cm	7.5cm	.038" (0.97mm)	Braided 145 PIG	5		5.5cm
<b>7779-Q3</b>	6	.051"/1.29mm	110cm	8.0cm	.038" (0.97mm)	Mod. Braided 145 PIG	10		5.5cm
<b>7779-T3</b>	6	.051"/1.29mm	110cm	7.5cm	.038" (0.97mm)	Braided 145 PIG	5		7.5cm
<b>7779-M3</b>	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided 145 PIG	5		5.5cm



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	#side ports	small pigtail curve	angle length
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	-------------	---------------------	--------------

**PERFORMA® 155 DEGREE CARDIAC PIGTAIL CATHETERS**

<b>7704-30</b>	4	.042"/1.07mm	110cm	6.5cm	.038" (0.97mm)	Braided 155 PIG	5		5.5cm
<b>7602-33</b>	5	.052"/1.32mm	110cm	7.5cm	.038" (0.97mm)	Non-Braided 155 PIG	5		7.5cm
<b>7529-33</b>	5	.046"/1.17mm	110cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		7.5cm
<b>7529-63</b>	5	.046"/1.17mm	125cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		7.5cm
<b>7776-D3</b>	5	.046"/1.17mm	110cm	4.6cm	.038" (0.97mm)	UHF 155 PIG	5	yes	5.5cm
<b>7776-43</b>	5	.046"/1.17mm	110cm	5.0cm	.038" (0.97mm)	Braided 155 PIG	5	yes	5.5cm
<b>7569-33</b>	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided 155 PIG	5		5.5cm
<b>7509-33</b>	6	.054"/1.37mm	110cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		5.5cm
<b>7509-63</b>	6	.054"/1.37mm	125cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		5.5cm
<b>7569-63</b>	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided 155 PIG	5		5.5cm
<b>7777-D3</b>	6	.059"/1.49mm	110cm	4.8cm	.038" (0.97mm)	Non-Braided 155 PIG	5	yes	5.5cm
<b>7777-43</b>	6	.054"/1.37mm	110cm	4.6cm	.038" (0.97mm)	Braided 155 PIG	5	yes	5.5cm
<b>7669-33</b>	7	.064"/1.62mm	110cm	6.5cm	.038" (0.97mm)	Braided 155 PIG	5		7.5cm
<b>7669-63</b>	7	.064"/1.62mm	125cm	6.5cm	.038" (0.97mm)	Braided 155 PIG	5		7.5cm

**SOFTOUCH® 155 DEGREE CARDIAC PIGTAIL CATHETERS**

<b>7779-73</b>	5	.052"/1.32mm	110cm	4.7cm	.038" (0.97mm)	Non-Braided 155 PIG	5		5.5cm
<b>7779-C3</b>	6	.051"/1.29mm	110cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		5.5cm
<b>7779-R3</b>	6	.051"/1.29mm	110cm	8.0cm	.038" (0.97mm)	Mod. Braided 155 PIG	10		5.5cm
<b>7779-U3</b>	6	.051"/1.29mm	110cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		7.5cm
<b>7779-N3</b>	6	.059"/1.49mm	110cm	8.0cm	.038" (0.97mm)	Non-Braided 155 PIG	5		5.5cm
<b>7780-33</b>	6	.051"/1.29mm	125cm	7.5cm	.038" (0.97mm)	Braided 155 PIG	5		7.5cm
<b>7222-33</b>	6	.059"/1.49mm	125cm	8.0cm	.038" (0.97mm)	Non-Braided 155 PIG	5		5.5cm



catalog number	French size	internal diameter	length	recommended guide wire	#side ports	pediatric adult
----------------	-------------	-------------------	--------	------------------------	-------------	-----------------

**PEDIATRIC / ADULT PIGTAIL CATHETERS - (FEP - TEFLON)\***

\* Teflon®-registered trademark of DuPont

<b>1310-21-0043</b>	3	.027"/.69mm	40cm	.021" (0.53mm)	4	P
<b>1310-21-0053</b>	3	.027"/.69mm	50cm	.021" (0.53mm)	4	P
<b>1310-21-0054</b>	4	.040"/1.02mm	50cm	.021" (0.53mm)	4	P
<b>1310-21-0064</b>	4	.040"/1.02mm	65cm	.021" (0.53mm)	4	P
<b>1310-25-0054</b>	4	.040"/1.02mm	50cm	.025" (0.64mm)	4	P
<b>1310-25-0064</b>	4	.040"/1.02mm	65cm	.025" (0.64mm)	4	P
<b>1310-25-0084</b>	4	.040"/1.02mm	80cm	.025" (0.64mm)	4	P
<b>1310-25-0065</b>	5	.046"/1.17mm	60cm	.025" (0.64mm)	4	P
<b>1310-25-0085</b>	5	.046"/1.17mm	80cm	.025" (0.64mm)	4	P
<b>1310-35-0086</b>	6	.056"/1.42mm	80cm	.035" (0.89mm)	4	P
<b>1308-25-0054</b>	4	.040"/1.02mm	50cm	.025" (0.64mm)	6	A
<b>1308-25-0085</b>	5	.046"/1.17mm	80cm	.025" (0.64mm)	6	A
<b>1308-35-0065</b>	5	.046"/1.17mm	65cm	.035" (0.89mm)	6	A
<b>1308-35-0085</b>	5	.046"/1.17mm	80cm	.035" (0.89mm)	6	A
<b>1308-35-0005</b>	5	.046"/1.17mm	100cm	.035" (0.89mm)	6	A
<b>1308-35-0086</b>	6	.056"/1.42mm	80cm	.035" (0.89mm)	6	A
<b>1308-35-0006</b>	6	.056"/1.42mm	100cm	.035" (0.89mm)	6	A
<b>1308-38-0006</b>	6	.056"/1.42mm	100cm	.038" (0.97mm)	6	A
<b>1308-35-0007</b>	7	.065"/1.65mm	100cm	.035" (0.89mm)	6	A
<b>1308-38-0007</b>	7	.065"/1.65mm	100cm	.038" (0.97mm)	6	A

catalog number	French size	internal diameter	length	length of tip	recommended guide wire	# side ports	description
----------------	-------------	-------------------	--------	---------------	------------------------	--------------	-------------

**SOFTOUCH® STANDARD PEDIATRIC AND INFANT PIGTAIL CATHETERS - (NYLON)**



<b>4406-11</b>	4	.041"/1.04mm	65cm	6.0cm	.035" (0.89mm)	5	standard non-braided
<b>4406-21</b>	4	.041"/1.04mm	80cm	6.0cm	.035" (0.89mm)	5	standard non-braided
<b>4408-A1</b>	4	.041"/1.04mm	65cm	5.5cm	.035" (0.89mm)	4	infant non-braided
<b>4408-B1</b>	4	.041"/1.04mm	80cm	5.5cm	.035" (0.89mm)	4	infant non-braided



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	# side ports
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--------------

**PERFORMA® PEDIATRIC JUDKINS RIGHT AND LEFT CATHETERS**

<b>7700-A0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 1.5	- bumper tip
<b>7700-B0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 2.0	- bumper tip
<b>7700-C0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 2.5	- bumper tip
<b>7700-D0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JL 3.0	- bumper tip
<b>7701-A0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 1.5	- bumper tip
<b>7701-B0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 2.0	- bumper tip
<b>7701-C0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 2.5	- bumper tip
<b>7701-D0</b>	4	.042"/1.07mm	70cm	2.0cm	.038" (0.97mm)	JR 3.0	- bumper tip
<b>7521-L3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 1.5	- bumper tip
<b>7521-M3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 2.0	- bumper tip
<b>7521-N3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 2.5	- bumper tip
<b>7521-P3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JL 3.0	- bumper tip
<b>7523-L3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 1.5	- bumper tip
<b>7523-M3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 2.0	- bumper tip
<b>7523-N3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 2.5	- bumper tip
<b>7523-P3</b>	5	.046"/1.17mm	70cm	2.0cm	.038" (0.97mm)	JR 3.0	- bumper tip



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	#side ports
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	-------------

### PERFORMA® MULTIPURPOSE CATHETERS

<b>7708-30</b>	4	.042"/1.07mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2 bumper tip
<b>7708-40</b>	4	.042"/1.07mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2 bumper tip
<b>7708-10</b>	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	- bumper tip
<b>7708-20</b>	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2 bumper tip
<b>7708-50</b>	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	- bumper tip
<b>7708-60</b>	4	.042"/1.07mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2 bumper tip
<b>7533-53</b>	5	.046"/1.17mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7533-63</b>	5	.046"/1.17mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7533-13</b>	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-
<b>7533-23</b>	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7533-33</b>	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-
<b>7533-43</b>	5	.046"/1.17mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7513-A3</b>	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI A1	-
<b>7513-13</b>	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7513-D3</b>	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI B1	-
<b>7513-33</b>	6	.054"/1.37mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7513-B3</b>	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-
<b>7513-23</b>	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7513-E3</b>	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-
<b>7513-43</b>	6	.054"/1.37mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7513-C3</b>	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI A1	-
<b>7513-53</b>	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7513-F3</b>	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI B1	-
<b>7513-63</b>	6	.054"/1.37mm	125cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7664-63</b>	7	.064"/1.62mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7664-83</b>	7	.064"/1.62mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7664-13</b>	7	.064"/1.62mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-
<b>7664-23</b>	7	.064"/1.62mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7664-33</b>	7	.064"/1.62mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-
<b>7664-43</b>	7	.064"/1.62mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	#side ports
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	-------------

**SOFTOUCH® MULTIPURPOSE CATHETERS**

<b>7407-23</b>	5	.044"/1.12mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7407-B3</b>	5	.044"/1.12mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7407-33</b>	5	.044"/1.12mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7407-C3</b>	5	.044"/1.12mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7407-M3</b>	5	.044"/1.12mm	100cm	5.8cm	.038" (0.97mm)	MULTI A1	-
<b>7407-R3</b>	5	.044"/1.12mm	100cm	5.8cm	.038" (0.97mm)	MULTI B1	-
<b>7407-63</b>	6	.051"/1.29mm	80cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7407-F3</b>	6	.051"/1.29mm	80cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7407-73</b>	6	.051"/1.29mm	100cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7407-G3</b>	6	.051"/1.29mm	100cm	5.8cm	.038" (0.97mm)	MULTI B2	2
<b>7407-83</b>	6	.051"/1.29mm	125cm	5.8cm	.038" (0.97mm)	MULTI A2	2
<b>7407-H3</b>	6	.051"/1.29mm	125cm	5.8cm	.038" (0.97mm)	MULTI B2	2



**SOFTOUCH® SONES CATHETERS**

<b>7481-13</b>	6	.051"/1.29mm	80cm	7.3cm	.038" (0.97mm)	SONES	4
<b>7481-23</b>	6	.051"/1.29mm	100cm	7.3cm	.038" (0.97mm)	SONES	4

**PERFORMA® SONES CATHETERS**

<b>7705-10</b>	4	.042"/1.07mm	80cm	7.3cm	.038" (0.97mm)	SONES	4	bumper tip
<b>7705-20</b>	4	.042"/1.07mm	100cm	7.3cm	.038" (0.97mm)	SONES	4	bumper tip
<b>7530-13</b>	5	.046"/1.17mm	80cm	7.3cm	.038" (0.97mm)	SONES	4	
<b>7530-23</b>	5	.046"/1.17mm	100cm	7.3cm	.038" (0.97mm)	SONES	4	
<b>7667-13</b>	7	.064"/1.62mm	80cm	7.5cm	.038" (0.97mm)	SONES	4	
<b>7667-23</b>	7	.064"/1.62mm	100cm	7.5cm	.038" (0.97mm)	SONES	4	



catalog number	French size	internal diameter	length	length of tip	recommended guide wire	description	
----------------	-------------	-------------------	--------	---------------	------------------------	-------------	--


**PERFORMA® INTERNAL MAMMARY LEFT - BRACHIAL / RADIAL CATHETERS**

<b>7242-10</b>	4	.042"/1.07mm	100cm	4.2cm	.038" (0.97mm)	IMLBR1	bumper tip
<b>7242-20</b>	4	.042"/1.07mm	100cm	4.2cm	.038" (0.97mm)	IMLBR2	bumper tip



**PERFORMA® LEFT CORONARY - BRACHIAL / RADIAL CATHETERS**

<b>7240-10</b>	4	.042"/1.07mm	100cm	3.4cm	.038" (0.97mm)	LBR1	bumper tip
<b>7240-20</b>	4	.042"/1.07mm	100cm	3.4cm	.038" (0.97mm)	LBR2	bumper tip
<b>7240-30</b>	4	.042"/1.07mm	100cm	6.0cm	.038" (0.97mm)	LBR3	bumper tip



**PERFORMA® RIGHT CORONARY - BRACHIAL / RADIAL CATHETER**

<b>7241-10</b>	4	.042"/1.07mm	100cm	3.4cm	.038" (0.97mm)	RBR	bumper tip
----------------	---	--------------	-------	-------	----------------	-----	------------



**PERFORMA® AMPLATZ LEFT CORONARY - BRACHIAL / RADIAL CATHETER**

<b>7703-B0</b>	4	.042"/1.07mm	100cm	1.5cm	.038" (0.97mm)	AL2	bumper tip
----------------	---	--------------	-------	-------	----------------	-----	------------




catalog number	French size	description	catheters and individual cat. no.
----------------	-------------	-------------	-----------------------------------

### PERFORMA° CARDIAC MULTIPACK

<b>7751-10</b>	4	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7704-10, 7700-20, 7701-20)	bumper tip
<b>7751-20</b>	4	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7704-20, 7700-20, 7701-20)	bumper tip
<b>7560-13</b>	5	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7529-11, 7521-21, 7523-21)	
<b>7560-A3</b>	5	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, PRC (7529-11, 7521-21, 7532-13)	
<b>7560-23</b>	5	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7529-23, 7521-21, 7523-21)	
<b>7560-C3</b>	5	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, PRC (7529-23, 7521-21, 7532-13)	
<b>7560-33</b>	5	CARDIAC MULTIPACK III	155° PIG (Braided), JL4.0, JR4.0 (7529-33, 7521-21, 7523-21)	
<b>7560-43</b>	5	CARDIAC MULTIPACK IV	PIG (Non-Braided), JL4.0, JR4.0 (7602-11, 7521-21, 7523-21)	
<b>7560-53</b>	5	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7602-23, 7521-21, 7523-21)	
<b>7660-13</b>	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7509-11, 7501-21, 7503-21)	
<b>7660-A3</b>	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, PRC (7509-11, 7501-21, 7512-13)	
<b>7660-23</b>	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7509-23, 7501-21, 7503-21)	
<b>7660-B3</b>	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, PRC (7509-23, 7501-21, 7512-13)	
<b>7660-33</b>	6	CARDIAC MULTIPACK III	155° PIG (Braided), JL4.0, JR4.0 (7509-33, 7501-21, 7503-21)	
<b>7660-43</b>	6	CARDIAC MULTIPACK IV	PIG (Non-Braided), JL4.0, JR4.0 (7569-11, 7501-21, 7503-21)	
<b>7660-53</b>	6	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7569-23, 7501-21, 7503-21)	
<b>7660-63</b>	6	CARDIAC MULTIPACK VI	155° PIG (Non-Braided), JL4.0, JR4.0 (7569-33, 7501-21, 7503-21)	

### SOFTOUCH° CARDIAC MULTIPACK

<b>7460-13</b>	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7779-A1, 7402-61, 7401-B1)
<b>7460-23</b>	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7779-B3, 7402-61, 7401-B1)
<b>7460-83</b>	6	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7779-M3, 7402-61, 7401-B1)




CARDIAC MULTIPACK

catalog number	French size	description	pigtail 110cm - JL & JR 100cm catheters	sheath	.035" (0.89mm) guide wire
----------------	-------------	-------------	---	--------	---------------------------

### CARDIAC MULTIPACK PLUS

<b>5MPSW-S</b>	5	CATHETERS, SHEATH AND GUIDE WIRE	PIG (Braided), JL4.0, JR4.0	5F PRELUDE*	150cm - J3
<b>5MPSW-145</b>	5	CATHETERS, SHEATH AND GUIDE WIRE	145° PIG ( Braided), JL4.0, JR4.0	5F PRELUDE*	150cm - J3
<b>5MPS-S</b>	5	CATHETERS, SHEATH	PIG (Braided), JL4.0, JR4.0	5F PRELUDE*	N/A
<b>5MPS-145</b>	5	CATHETERS, SHEATH	145° PIG ( Braided), JL4.0, JR4.0	5F PRELUDE*	N/A
<b>6MPSW-S</b>	6	CATHETERS, SHEATH AND GUIDE WIRE	PIG (Braided), JL4.0, JR4.0	6F PRELUDE*	150cm - J3
<b>6MPSW-145</b>	6	CATHETERS, SHEATH AND GUIDE WIRE	145° PIG ( Braided), JL4.0, JR4.0	6F PRELUDE*	150cm - J3
<b>6MPS-S</b>	6	CATHETERS, SHEATH	PIG (Braided), JL4.0, JR4.0	6F PRELUDE*	N/A
<b>6MPS-145</b>	6	CATHETERS, SHEATH	145° PIG ( Braided), JL4.0, JR4.0	6F PRELUDE*	N/A



**CARDIAC MULTIPACK PLUS**

\* Sheath, Dilator, 50cm .035" (0.97mm) and J3 Guide Wire