

Guide Catheters

Coronary angioplasty requires the use of a guiding catheter to place the balloon within the vasculature. Once in place, the guiding catheter acts as a conduit for wires and catheters, and radiopaque dye that is used to provide fluoroscopic visualization during the procedure.



Guide Catheters

GUIDE CATHETERS





GUIDE CATHETERS

D

TRAX®

Judkins Left	D.1
Judkins Left Short Tip	D.2
Judkins Left Flexible Tip	D.2
Judkins Left Flexible Short Tip	D.3
Left Support	D.3
Judkins Right	D.4
Judkins Right Short Tip	D.4
Judkins Right Flexible Tip	D.5
Judkins Right Flexible Short Tip	D.5
Right Support	D.6
Amplatz Left	D.6
Amplatz Left Short Tip	D.7
Amplatz Right	D.7, 8
Amplatz Right Short Tip	D.8
Progressive Right	D.9
Progressive Right II	D.9
Progressive Left I	D.9
Coronary Bypass	D.10
Internal Mammary Artery	D.10
Multipurpose	D.11
Hockey Stick	D.12
Renal Curves	D.13
Cobra Curves	D.13

catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------




TRAX™ JUDKINS LEFT

TRAX 7 FRENCH

7821-A0	Judkins Left 3.5	.074"/1.88mm	.092"/2.34mm	100cm	-
7821-E0	Judkins Left 3.5 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7821-B0	Judkins Left 4.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7821-F0	Judkins Left 4.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7821-J0	Judkins Left 4.5	.074"/1.88mm	.092"/2.34mm	100cm	-
7821-K0	Judkins Left 4.5 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7821-C0	Judkins Left 5.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7821-G0	Judkins Left 5.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7821-D0	Judkins Left 6.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7821-H0	Judkins Left 6.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7831-L0	Judkins Left 3.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-P0	Judkins Left 3.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-A0	Judkins Left 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-F0	Judkins Left 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-B0	Judkins Left 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-G0	Judkins Left 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-C0	Judkins Left 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-H0	Judkins Left 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-D0	Judkins Left 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-J0	Judkins Left 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-E0	Judkins Left 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-K0	Judkins Left 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ JUDKINS LEFT SHORT TIPS




TRAX 8 FRENCH

7831-10	Judkins Left 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-50	Judkins Left 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-20	Judkins Left 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-60	Judkins Left 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-M0	Judkins Left 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-N0	Judkins Left 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-30	Judkins Left 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-70	Judkins Left 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7831-40	Judkins Left 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7831-80	Judkins Left 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2

TRAX™ JUDKINS LEFT FLEXIBLE TIPS

TRAX 8 FRENCH

7854-10	Judkins Left 3.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7854-A0	Judkins Left 3.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7854-20	Judkins Left 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7854-B0	Judkins Left 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7854-30	Judkins Left 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7854-C0	Judkins Left 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7854-40	Judkins Left 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7854-D0	Judkins Left 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7854-50	Judkins Left 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7854-E0	Judkins Left 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7854-60	Judkins Left 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7854-F0	Judkins Left 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ JUDKINS LEFT FLEXIBLE SHORT TIPS

TRAX 8 FRENCH

7862-10	Judkins Left 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7862-A0	Judkins Left 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7862-20	Judkins Left 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7862-B0	Judkins Left 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7862-30	Judkins Left 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7862-C0	Judkins Left 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7862-40	Judkins Left 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7862-D0	Judkins Left 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7862-50	Judkins Left 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7862-E0	Judkins Left 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2




TRAX™ LEFT SUPPORT

TRAX 7 FRENCH

7858-10	Left Support 3.5	.074"/1.88mm	.092"/2.34mm	100cm	-
7858-A0	Left Support 3.5 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7858-20	Left Support 4.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7858-B0	Left Support 4.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7858-30	Left Support 4.5	.074"/1.88mm	.092"/2.34mm	100cm	-
7858-C0	Left Support 4.5 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7858-40	Left Support 5.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7858-D0	Left Support 5.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7856-10	Left Support 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7856-A0	Left Support 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7856-20	Left Support 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7856-B0	Left Support 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7856-40	Left Support 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7856-D0	Left Support 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7856-30	Left Support 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7856-C0	Left Support 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ JUDKINS RIGHT

TRAX 7 FRENCH

7820-AO	Judkins Right 3.5	.074"/1.88mm	.092"/2.34mm	100cm	-
7820-EO	Judkins Right 3.5 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7820-BO	Judkins Right 4.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7820-FO	Judkins Right 4.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7820-JO	Judkins Right 4.5	.074"/1.88mm	.092"/2.34mm	100cm	-
7820-KO	Judkins Right 4.5 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7820-CO	Judkins Right 5.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7820-GO	Judkins Right 5.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7820-DO	Judkins Right 6.0	.074"/1.88mm	.092"/2.34mm	100cm	-
7820-HO	Judkins Right 6.0 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7830-LO	Judkins Right 3.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-RO	Judkins Right 3.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-AO	Judkins Right 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-FO	Judkins Right 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-BO	Judkins Right 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-MO	Judkins Right 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	90cm	2
7830-GO	Judkins Right 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-CO	Judkins Right 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-HO	Judkins Right 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-DO	Judkins Right 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-JO	Judkins Right 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-EO	Judkins Right 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-KO	Judkins Right 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



TRAX™ JUDKINS RIGHT SHORT TIPS

TRAX 8 FRENCH

7830-10	Judkins Right 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-50	Judkins Right 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-20	Judkins Right 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-60	Judkins Right 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-N0	Judkins Right 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-P0	Judkins Right 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-30	Judkins Right 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-70	Judkins Right 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7830-40	Judkins Right 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7830-80	Judkins Right 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2

catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ JUDKINS RIGHT FLEXIBLE TIPS


TRAX 8 FRENCH

7853-10	Judkins Right 3.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7853-A0	Judkins Right 3.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7853-20	Judkins Right 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7853-B0	Judkins Right 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7853-30	Judkins Right 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7853-C0	Judkins Right 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7853-40	Judkins Right 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7853-D0	Judkins Right 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7853-50	Judkins Right 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7853-E0	Judkins Right 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7853-60	Judkins Right 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7853-F0	Judkins Right 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2

TRAX™ JUDKINS RIGHT FLEXIBLE SHORT TIPS

TRAX 8 FRENCH

7861-10	Judkins Right 3.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7861-A0	Judkins Right 3.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7861-20	Judkins Right 4.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7861-B0	Judkins Right 4.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7861-30	Judkins Right 4.5	.089"/2.26mm	.104"/2.64mm	100cm	-
7861-C0	Judkins Right 4.5 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7861-40	Judkins Right 5.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7861-D0	Judkins Right 5.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7861-50	Judkins Right 6.0	.089"/2.26mm	.104"/2.64mm	100cm	-
7861-E0	Judkins Right 6.0 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ RIGHT SUPPORT

TRAX 7 FRENCH

7857-10	Right Support 1	.074"/1.88mm	.092"/2.34mm	100cm	-
7857-A0	Right Support 1 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7857-20	Right Support 2	.074"/1.88mm	.092"/2.34mm	100cm	-
7857-B0	Right Support 2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7857-30	Right Support 3	.074"/1.88mm	.092"/2.34mm	100cm	-
7857-C0	Right Support 3 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7855-10	Right Support 1	.089"/2.26mm	.104"/2.64mm	100cm	-
7855-A0	Right Support 1 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7855-20	Right Support 2	.089"/2.26mm	.104"/2.64mm	100cm	-
7855-B0	Right Support 2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7855-30	Right Support 3	.089"/2.26mm	.104"/2.64mm	100cm	-
7855-C0	Right Support 3 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



TRAX™ AMPLATZ LEFT

TRAX 7 FRENCH

7819-50	Amplatz Left 1	.074"/1.88mm	.092"/2.34mm	100cm	-
7819-A0	Amplatz Left 1 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7819-60	Amplatz Left 2	.074"/1.88mm	.092"/2.34mm	100cm	-
7819-B0	Amplatz Left 2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7819-70	Amplatz Left 3	.074"/1.88mm	.092"/2.34mm	100cm	-
7819-C0	Amplatz Left 3 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7819-80	Amplatz Left 4	.074"/1.88mm	.092"/2.34mm	100cm	-
7819-D0	Amplatz Left 4 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7833-10	Amplatz Left 1	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-A0	Amplatz Left 1 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7833-20	Amplatz Left 2	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-B0	Amplatz Left 2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7833-30	Amplatz Left 3	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-C0	Amplatz Left 3 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7833-40	Amplatz Left 4	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-D0	Amplatz Left 4 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ AMPLATZ LEFT SHORT TIPS

TRAX 7 FRENCH

7864-10	Amplatz Left 1	.074"/1.88mm	.092"/2.34mm	100cm	-
7864-A0	Amplatz Left 1 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7864-20	Amplatz Left 2	.074"/1.88mm	.092"/2.34mm	100cm	-
7864-B0	Amplatz Left 2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7864-30	Amplatz Left 3	.074"/1.88mm	.092"/2.34mm	100cm	-
7864-C0	Amplatz Left 3 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7864-40	Amplatz Left 4	.074"/1.88mm	.092"/2.34mm	100cm	-
7864-D0	Amplatz Left 4 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2


TRAX 8 FRENCH

7833-50	Amplatz Left 1	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-E0	Amplatz Left 1 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7833-60	Amplatz Left 2	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-F0	Amplatz Left 2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7833-70	Amplatz Left 3	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-G0	Amplatz Left 3 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7833-80	Amplatz Left 4	.089"/2.26mm	.104"/2.64mm	100cm	-
7833-H0	Amplatz Left 4 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2

TRAX™ AMPLATZ RIGHT

TRAX 7 FRENCH

7814-40	Amplatz Right 1	.074"/1.88mm	.092"/2.34mm	100cm	-
7814-A0	Amplatz Right 1 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7814-50	Amplatz Right 2	.074"/1.88mm	.092"/2.34mm	100cm	-
7814-B0	Amplatz Right 2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7814-60	Amplatz Right 3	.074"/1.88mm	.092"/2.34mm	100cm	-
7814-C0	Amplatz Right 3 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7814-70	Amplatz Right Modified	.074"/1.88mm	.092"/2.34mm	100cm	-
7814-H0	Amplatz Right Modified with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ AMPLATZ RIGHT

TRAX 8 FRENCH

7832-10	Amplatz Right 1	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-A0	Amplatz Right 1 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7832-20	Amplatz Right 2	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-B0	Amplatz Right 2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7832-30	Amplatz Right 3	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-C0	Amplatz Right 3 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7832-70	Amplatz Right Modified	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-G0	Amplatz Right Modified with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2

TRAX™ AMPLATZ RIGHT SHORT-TIPS

TRAX 7 FRENCH

7863-10	Amplatz Right 1	.074"/1.88mm	.092"/2.34mm	100cm	-
7863-A0	Amplatz Right 1 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7863-20	Amplatz Right 2	.074"/1.88mm	.092"/2.34mm	100cm	-
7863-B0	Amplatz Right 2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7863-30	Amplatz Right 3	.074"/1.88mm	.092"/2.34mm	100cm	-
7863-C0	Amplatz Right 3 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7863-40	Amplatz Right Modified	.074"/1.88mm	.092"/2.34mm	100cm	-
7863-D0	Amplatz Right Modified with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7832-40	Amplatz Right 1	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-D0	Amplatz Right 1 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7832-50	Amplatz Right 2	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-E0	Amplatz Right 2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7832-60	Amplatz Right 3	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-F0	Amplatz Right 3 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7832-80	Amplatz Right Modified	.089"/2.26mm	.104"/2.64mm	100cm	-
7832-H0	Amplatz Right Modified with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2

catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------


TRAX® PROGRESSIVE RIGHT

TRAX 7 FRENCH

7828-10	Progressive Right	.074"/1.88mm	.092"/2.34mm	100cm	-
7828-A0	Progressive Right with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7829-20	Progressive Right	.089"/2.26mm	.104"/2.64mm	100cm	-
7829-C0	Progressive Right with side ports	.089"/2.26mm	.104"/2.64mm	90cm	2
7829-B0	Progressive Right with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2




TRAX® PROGRESSIVE RIGHT II

TRAX 7 FRENCH

7809-30	Progressive Right II	.074"/1.88mm	.092"/2.34mm	100cm	-
7809-C0	Progressive Right II with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7809-10	Progressive Right II	.089"/2.26mm	.104"/2.64mm	100cm	-
7809-A0	Progressive Right II with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2




TRAX® PROGRESSIVE LEFT I

TRAX 7 FRENCH

7808-30	Progressive Left I	.074"/1.88mm	.092"/2.34mm	100cm	-
7808-C0	Progressive Left I with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7808-10	Progressive Left I	.089"/2.26mm	.104"/2.64mm	100cm	-
7808-A0	Progressive Left I with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ CORONARY BYPASS

TRAX 7 FRENCH

7807-50	Right Coronary Bypass	.074"/1.88mm	.092"/2.34mm	55cm	-
7807-30	Right Coronary Bypass	.074"/1.88mm	.092"/2.34mm	100cm	-
7807-40	Left Coronary Bypass	.074"/1.88mm	.092"/2.34mm	100cm	-

TRAX 8 FRENCH

7807-C0	Right Coronary Bypass	.089"/2.26mm	.104"/2.64mm	55cm	-
7807-A0	Right Coronary Bypass	.089"/2.26mm	.104"/2.64mm	100cm	-
7807-B0	Left Coronary Bypass	.089"/2.26mm	.104"/2.64mm	100cm	-



TRAX™ INTERNAL MAMMARY ARTERY

TRAX 7 FRENCH

7463-20	Internal Mammary Artery	.074"/1.88mm	.092"/2.34mm	100cm	-
----------------	-------------------------	--------------	--------------	-------	---

TRAX 8 FRENCH

7463-30	Internal Mammary Artery	.089"/2.26mm	.104"/2.64mm	100cm	-
----------------	-------------------------	--------------	--------------	-------	---



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------


TRAX™ MULTIPURPOSE

TRAX 7 FRENCH

7860-E0	Multipurpose A1	.074"/1.88mm	.092"/2.34mm	55cm	-
7860-G0	Multipurpose A2 with side ports	.074"/1.88mm	.092"/2.34mm	55cm	2
7860-F0	Multipurpose B1	.074"/1.88mm	.092"/2.34mm	55cm	-
7860-H0	Multipurpose B2 with side ports	.074"/1.88mm	.092"/2.34mm	55cm	2
7860-50	Multipurpose A1	.074"/1.88mm	.092"/2.34mm	65cm	-
7860-70	Multipurpose A2 with side ports	.074"/1.88mm	.092"/2.34mm	65cm	2
7860-60	Multipurpose B1	.074"/1.88mm	.092"/2.34mm	65cm	-
7860-80	Multipurpose B2 with side ports	.074"/1.88mm	.092"/2.34mm	65cm	2
7815-50	Multipurpose A1	.074"/1.88mm	.092"/2.34mm	100cm	-
7815-70	Multipurpose A2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2
7815-60	Multipurpose B1	.074"/1.88mm	.092"/2.34mm	100cm	-
7815-80	Multipurpose B2 with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7815-J0	Multipurpose A1	.089"/2.26mm	.104"/2.64mm	55cm	-
7815-L0	Multipurpose A2 with side ports	.089"/2.26mm	.104"/2.64mm	55cm	2
7815-K0	Multipurpose B1	.089"/2.26mm	.104"/2.64mm	55cm	-
7815-M0	Multipurpose B2 with side ports	.089"/2.26mm	.104"/2.64mm	55cm	2
7815-E0	Multipurpose A1	.089"/2.26mm	.104"/2.64mm	65cm	-
7815-G0	Multipurpose A2 with side ports	.089"/2.26mm	.104"/2.64mm	65cm	2
7815-F0	Multipurpose B1	.089"/2.26mm	.104"/2.64mm	65cm	-
7815-H0	Multipurpose B2 with side ports	.089"/2.26mm	.104"/2.64mm	65cm	2
7815-A0	Multipurpose A1	.089"/2.26mm	.104"/2.64mm	100cm	-
7815-C0	Multipurpose A2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2
7815-B0	Multipurpose B1	.089"/2.26mm	.104"/2.64mm	100cm	-
7815-D0	Multipurpose B2 with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------

TRAX™ HOCKEY STICK

TRAX 7 FRENCH

7845-JO	Hockey Stick I	.074"/1.88mm	.092"/2.34mm	55cm	-
7845-KO	Hockey Stick I with side ports	.074"/1.88mm	.092"/2.34mm	55cm	2
7845-CO	Hockey Stick I	.074"/1.88mm	.092"/2.34mm	65cm	-
7845-DO	Hockey Stick I with side ports	.074"/1.88mm	.092"/2.34mm	65cm	2
7845-AO	Hockey Stick I	.074"/1.88mm	.092"/2.34mm	100cm	-
7845-BO	Hockey Stick I with side ports	.074"/1.88mm	.092"/2.34mm	100cm	2

TRAX 8 FRENCH

7845-FO	Hockey Stick I	.089"/2.26mm	.104"/2.64mm	45cm	-
7845-EO	Hockey Stick I with side ports	.089"/2.26mm	.104"/2.64mm	45cm	2
7845-GO	Hockey Stick I	.089"/2.26mm	.104"/2.64mm	55cm	-
7845-HO	Hockey Stick I with side ports	.089"/2.26mm	.104"/2.64mm	55cm	2
7845-7O	Hockey Stick I	.089"/2.26mm	.104"/2.64mm	65cm	-
7845-9O	Hockey Stick I with side ports	.089"/2.26mm	.104"/2.64mm	65cm	2
7845-5O	Hockey Stick I	.089"/2.26mm	.104"/2.64mm	100cm	-
7845-3O	Hockey Stick I with side ports	.089"/2.26mm	.104"/2.64mm	100cm	2



catalog number	description	inner diameter	outer diameter	length	# of side ports
----------------	-------------	----------------	----------------	--------	-----------------




TRAX™ RENAL CURVES

TRAX 7 FRENCH

7859-30	Renal I	.074"/1.88mm	.092"/2.34mm	55cm	-
7859-60	Renal II	.074"/1.88mm	.092"/2.34mm	55cm	-
7859-C0	Renal Double Curve	.074"/1.88mm	.092"/2.34mm	55cm	-
7834-A0	Renal I	.074"/1.88mm	.092"/2.34mm	65cm	-
7805-10	Renal II	.074"/1.88mm	.092"/2.34mm	65cm	-
7806-10	Renal Double Curve	.074"/1.88mm	.092"/2.34mm	65cm	-

TRAX 8 FRENCH

7859-10	Renal I	.089"/2.26 mm	.104"/2.64mm	55cm	-
7859-40	Renal II	.089"/2.26 mm	.104"/2.64mm	55cm	-
7859-A0	Renal Double Curve	.089"/2.26 mm	.104"/2.64mm	55cm	-
7834-10	Renal I	.089"/2.26 mm	.104"/2.64mm	65cm	-
7805-40	Renal II	.089"/2.26 mm	.104"/2.64mm	65cm	-
7806-40	Renal Double Curve	.089"/2.26 mm	.104"/2.64mm	65cm	-



TRAX™ COBRA CURVES

TRAX 7 FRENCH

7822-10	Cobra 2	.074"/1.88mm	.092"/2.34mm	55cm	-
----------------	---------	--------------	--------------	------	---

TRAX 8 FRENCH

7823-10	Cobra 2	.089"/2.26mm	.104"/2.64mm	55cm	-
7823-A0	Cobra 2 with side ports	.089"/2.26mm	.104"/2.64mm	55cm	2

