

*Sheath Introducers are used to
the access through which guide
and catheters, are passed into the
vasculature and incorporate a
to minimize bleeding and a side
port for drug delivery, pressure monitoring
and standard flushing.*




Introducer Sheaths & Accessories

Introducers & Accessories

INTRODUCERS AND ACCESSORIES

Introducer Sheath & Accessories




INTRODUCERS & ACCESSORIES

G

PRELUDE™ Sheath Introducers and Accessories	G.1
MERIT® MAK Introducers and Accessories	G.2, 3
DIALEASE™ Introducer Sheaths	G.4
Vessel Dilators	G.4

PRELUDE™ SHEATH INTRODUCERS

The Prelude™ Sheath Introducer is packaged in a holster that keeps all components secure, especially the guide wire. Components can be flushed while in the holster. Guide wires are double ended, stainless steel, 0.97mm with a floppy 3mm J end and a floppy straight end. Product includes a vessel dilator, sheath introducer and guide wire (also available without guide wire).



catalog number	French size	length	guide wire diameter/length	hub color
SHEATH INTRODUCERS				
PSI-4F-11-035	4F	11cm	0.035" (0.89mm) x 50cm	red
PSI-5F-11-038	5F	11cm	0.038" (0.97mm) x 50cm	gray
PSI-6F-11-038	6F	11cm	0.038" (0.97mm) x 50cm	green
PSI-7F-11-038	7F	11cm	0.038" (0.97mm) x 50cm	orange
PSI-8F-11-038	8F	11cm	0.038" (0.97mm) x 50cm	blue
PSI-4F-11	4F	11cm	N/A	red
PSI-5F-11	5F	11cm	N/A	gray
PSI-6F-11	6F	11cm	N/A	green
PSI-7F-11	7F	11cm	N/A	orange
PSI-8F-11	8F	11cm	N/A	blue

Packaged in holster; 5 per box

catalog number	French size	length	recommended sheath size	hub color
OBTURATORS				
OBT-4F-11	4F	13cm	5F x 11cm	red
OBT-5F-11	5F	13cm	6F x 11cm	gray
OBT-6F-11	6F	13cm	7F x 11cm	green
OBT-7F-11	7F	13cm	8F x 11cm	orange



Packaged in pouch; 10 per box

catalog number	length	guide wire diameter
GUIDE WIRES		
SW35DE50SST	50cm	0.035" (0.89mm)
SW38DE50SST	50cm	0.038" (0.97mm)

Packaged in hoop; 10 per box

MERIT® MAK (Mini Access Kits)

The Merit MAK consists of a 4F or 5F coaxial introducer/dilator pair, 21 gauge introducer needle, and 0.018" (0.046mm) x 40cm guide wire. 21G needles are available echo-enhanced (4cm & 7cm) and non-echo enhanced (7cm). Needle and coaxial pair are pouched individually and contained in an outer kit pouch with the guide wire (guide wire is in a hoop). Both 4F and 5F products accept an 0.035" (0.89mm) or 0.038" (0.97mm) guide wire.





catalog number	French size	length	guide wire diameter/length	guide wire core/tip	needle gauge/length	qty/box
STANDARD 7cm NEEDLE KITS						
MAK401	4F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
MAK501	5F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
MAK401N	4F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
MAK501N	5F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
MAK405	4F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
MAK505	5F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
ECHO ENHANCED 7cm NEEDLE KITS						
MAK401-E	4F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
MAK501-E	5F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
MAK401N-E	4F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
MAK501N-E	5F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
MAK405-E	4F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
MAK505-E	5F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
ECHO ENHANCED 4cm NEEDLE KITS						
MAK401-PED	4F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	10
MAK501-PED	5F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	10
MAK401N-PED	4F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	10
MAK501N-PED	5F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	10
MAK405-PED	4F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 4cm	10
MAK505-PED	5F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 4cm	10

SS = Stainless Steel PT = Platinum Tip NT = Nitinol

MERIT® S-MAK (Stiffened Dilator Mini Access Kits)

The Merit S-MAK consists of a 4F or 5F coaxial introducer/stiffened dilator pair, 21 gauge echo enhanced introducer needle, and 0.018" (0.46mm) x 40cm guide wire. Needle and coaxial pair are pouched individually and contained in an outer kit pouch with the guide wire (guide wire is in a hoop). Both 4F and 5F products accept an 0.035" (0.89mm) or 0.038" (0.97mm) guide wire. Packaged 10 units per box.



catalog number	French size	length	guide wire diameter/length	guide wire core/tip	needle gauge/length	qty/box
STIFFENED DIALATOR KITS - ECHO ENHANCED 7cm NEEDLE						
S-MAK401	4F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
S-MAK501	5F	10cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 7cm	10
S-MAK401N	4F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
S-MAK501N	5F	10cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	10
S-MAK405	4F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10
S-MAK505	5F	10cm	0.018" (0.46mm) x 40cm	SS/PT	21G x 7cm	10

SS = Stainless Steel PT = Platinum Tip NT = Nitinol

GUIDE WIRES

The wires are contained in a hoop & individually pouched; 10 units per box.

Stainless steel guide wire with stainless steel tip

MAK001 0.018" (0.46mm) x 40cm

Nitinol guide wire with platinum tip

MAK001N 0.018" (0.46mm) x 40cm

MAK001N-60 0.018" (0.46mm) x 60cm

Stainless steel guide wire with platinum tip

MAK005 0.018" (0.46mm) x 40cm

MAK005-60 0.018" (0.46mm) x 60cm

NEEDLES


The needles are individually pouched; 25 units per box.

MAK21T71W	21G	7cm	
MAK21T71W-E	21G	7cm	echo enhanced
MAK21T41W-E	21G	4cm	echo enhanced

catalog number	French size	type	length	recommended guide wire	qty/box
----------------	-------------	------	--------	------------------------	---------

DIALEASE™ INTRODUCER SHEATHS

7011-A3 *	5	Arterial (Red)	4cm	.038" (0.97mm)	5
7011-13 *	5	Venous (Blue)	4cm	.038" (0.97mm)	5
7012-A3 *	6	Arterial (Red)	4cm	.038" (0.97mm)	5
7012-13 *	6	Venous (Blue)	4cm	.038" (0.97mm)	5
7013-A3 *	7	Arterial (Red)	4cm	.038" (0.97mm)	5
7013-13 *	7	Venous (Blue)	4cm	.038" (0.97mm)	5
7014-A3 *	8	Arterial (Red)	4cm	.038" (0.97mm)	5
7014-13 *	8	Venous (Blue)	4cm	.038" (0.97mm)	5



Radiopaque sheath with hemostasis valve and sideport with snap clamp and attached female luer fitting.

Radiopaque Vessel Dilator.

catalog number	French size	length	recommended guide wire	qty/box
----------------	-------------	--------	------------------------	---------

VESSEL DILATORS

3339-21	5	20cm	.035" (0.89mm)	10
3339-31	4	20cm	.035" (0.89mm)	10
3339-41	6	20cm	.038" (0.97mm)	10
3339-51	5	20cm	.038" (0.97mm)	10

* Not available in countries where CE mark is required.